

Neuropathy Analyser™ VIBROTHERM Dx®

(Combined Vibration Perception Threshold and Thermal Perception Threshold Detector)

An Instrument to detect small fibre and large fibre neuropathy..



Failure to perceive variances in temperature in the extremities is the most decisive early symptom of distal symmetric polyneuropathy, a typical ailment accompanying diabetes mellitus (Ziegler 1988).

Small-fiber Test provides an easy-to-use and scientifically validated measure of warm, cool and heat pain thermal sensory thresholds, all clinically useful determinants in the evaluation of diabetic neuropathy.

Vibration Perception Threshold (VPT) has been shown to be strongly associated with foot ulceration. VPT determination by using a Digital Biothesiometer has been used to identify peripheral sensory neuropathy and subjects at risk of foot. Digital Biothesiometer helps us to quantitate the threshold and monitor progressive changes or trends on following up testing.

Sensory neuropathy increases the risk of foot ulcerations by seven folds and peripheral arterial disease (PAD) by three folds in people with diabetes. Peripheral neuropathy is the major causal factors in the development of foot ulcerations among diabetic subjects.

Early detection helps better prevention. Prevention is better than no cure. 50% of all non-traumatic amputations occur in diabetes who is a high risk group.

FEATURES:

- ◆ Full solid state design with PC enabled
- ◆ Easy Tool to quantify Neuropathy
- ◆ Digital Biothesiometer and Thermometry HCP combined
- ◆ Digital 0 ~ 50 Volts Vibration output
- ◆ 1 ~ 50° C temperature output range
- ◆ 1° C rate of change per second
- ◆ Remote keypad helps easy operation
- ◆ Single probe delivers both Hot & Cold temperatures
- ◆ Automatic interpretation based on set values
- ◆ Early detection of Small Fibre Neuropathy
- ◆ Software supports Windows XP(SP3) / Win 7 / Win 8 / Win10 OS
- ◆ Patient report storage, export, pdf options available
- ◆ USB interface





Care for Life

Patient Name : Mohammed Sajan.O.A

Patient ID : 2

User ID : 1

G. D. Diabetes Institute and speciality clinics pvt. ltd

Age & Sex : 23 / Male

Date / Time : 28-11-2008 / 11:34:24 AM

Referred By: Dr.Sujay Majumdar

Right Foot

Left Foot

Right Foot

Left Foot

Mono Filament

Vibration







**Impression

Average: 6 - Normal

6 - Normal

Cold

Hot



Average: 22 - Normal

22.3 - Normal

Average: 39 - Normal

36 - Normal

**Impression

**Impression

Dr. Shaibal Chakraborty
Consultant Diabetologist

For more information contact



Diabetic Foot Care India Pvt Limited

No.18/1, Kannappanagar, 3rd Main Road, Thiruvannamipur, Chennai-41, INDIA.

Tel - 91-44-4356 4129 / 4359 0510 / 93806 21607

E-mail: mesmedi@gmail.com; elango@mesmedi.com

URL: www.diabeticfootcareindia.com; www.diabeticneuropathy.co.in